

This **Medicine Storage Box** is streamline, easy to carry and convenient for storing your medicine and medical supplies inside. Comes with a coded lock design to protect your goods and also safe for keeping chemicals & medicines out of reach from children.

Function

- 1 Medical Supplies Area
- 2 Bottled area
- 3 Code lock



Product Code: **TEB0353**

Precautions and Maintenance

1. To avoid damage to this box, avoid dropping from a great height.
2. Do not lift from handle if box is unlocked.
3. The initial password for the coded lock is 0000.
4. For safety considerations, please do not store valuable items inside.
5. Keep the product indoors.

Steps to change the coded lock





1. The initial code is set as 0000.

2. Switch the lever from point B to point A, after confirming the initial code.

3. If you wish to set your own pin code, please do so now.

4. Then switch the lever back from point A to point B to set the code.

5. Please remember your pin code, and if you wish to change again follow steps 2-4 as above.



Kind Reminder: For the safety and consideration of others, this box does not have any position for an outer unlock device in its structural design and it is not easy to dismantle after assembly. Please pay full attention to each step of code setting. Especially make sure the 4 digits lock code set is exactly in a parallel line before putting the lever in to position B. And do not forget your new set code to avoid any inconvenience.